#### Where's My Peak? Separating Truth from Fiction in Measurements of Nanoparticles Franklin Monzon, Lew Brown, Mac Bailey, Ngoc Do, Kaitlyn Huynh, Andrew Cleland, Peter Meinhold, Jean-Luc Fraikin www.NanoparticleAnalyzer.com Spectradyne LLC, 23875 Madison St. Suite A, Torrance, CA 90505 info@spectradynellc.com

#### Problem Statement

Regardless of an instrument's underlying principle of operation, all nanoparticle analysis instruments (and, indeed, instruments generally) are eventually limited by their intrinsic sensitivity and by instrument noise. This poster attempts to clarify relevant considerations in the detection of submicron particles near an instrument's limit of detection, providing a guide for researchers facing this kind of small-signal challenge.



#### Example 1: Extracellular Vesicle (EV) Measurement with the nCS1



## 3.0E+04 2.5E+04



#### Example 2: Protein Aggregation Measurement with the nCS1 and NTA

#### Sensitivity differences between instruments result in different detection thresholds, as is clear in these protein aggregation measurements.



#### Example 3: Exosome Quantification with the nCS1, NTA, and Cryo TEM

#### A three-way comparison of an exosome measurement illustrates the NTA false peak issue clearly. Cryo TEM and nCS1 agree closely.

NTA measurement of urinary exosomes pooled from five donors shows a pronounced peak near 140 nm. In this plot the y-axis has a linear scale.





MRPS is a state-of-the-art microfluidic implementation of Resistive Pulse Sensing, aka the Coulter Principle. In RPS, the electrical resistance of a conducting fluid is monitored as particles flow through a constriction (see illustration below), thereby blocking the flow of ions and temporarily increasing the resistance. Maximum resistance modulation is obtained when the particle is midway through the constriction. The size of the resistance spike is proportional to particle volume, regardless of particle material, and the transit time gives the fluid flow rate, which en-NC ables absolute concentration measurements. Flow rate

There is no dependence on the index of refraction of the particle and, because particles are detect-Time (µs) ed and sized individually, high resolution measurements are obtained. Spectradyne's proprietary MRPS technology utilizes disposable cartridges to greatly improve ease-of-use and reduce measurement time compared to other nanoparticle implementations of RPS.

#### Key Features of the nCS1 and MRPS:

### Recommendations

As with any measurement, researchers must understand instrumentation limitations prior to drawing conclusions about nanoparticle distributions. Detection thresholds, which may vary by particle type (for light scattering-based instruments, at least) need to be understood and any peaks in the distribution near the threshold should be greeted with skepticism until orthogonal methods can be employed. Failing that, dilution series can be run to confirm that detection events scale with anticipated concentration changes. In addition, checks in data analysis procedures should be in place to screen out false positives. Although taking care to confirm that peaks are real requires time and money, this investment is far less expensive than the consequences of drawing erroneous conclusions about one's formulations, which can lead to poor therapeutic and diagnostic outcomes.



→ nCS1 1(

NTA

# **Spectradyne**<sup>®</sup>



#### **Particle Analysis**

#### **Overview of Spectradyne's Technology**



• A truly orthogonal method to DLS, NTA, etc. • Sizing range: 50 nm - 10 μm diameter. • Peak sizing resolution of 3% or better. Absolute concentration measurements. Concentration range: 10<sup>5</sup> to 10<sup>12</sup> NPs/mL • All particle materials.

• Measurements unbiased by polydispersity. • Total sample analysis in minutes.



analysis using a disposable microfluidic cartridge (right), which pr ing requirement

**NSF & NIH SBIR** Phase II Awardee



NIH Biomedical In





